



CompactCath®

Catheterization Made Easy

Pre-Lubricated, Discreet, Touch Free.
100% DEHP, BPA, & Latex Free.

| Size | Tip | Item Number |
|-------|----------|-------------|
| 12 FR | Straight | 40788-1612 |
| 14 FR | Straight | 40788-1614 |
| 14 FR | Coude | 40788-1614C |



CompactCath® OneCath

Same High Quality At An Affordable Price Point

100% DEHP, BPA, and Latex Free

| Size | Tip | Item Number |
|-------|----------|-------------|
| 10 FR | Straight | 201-1610 |
| 12 FR | Straight | 201-1612 |
| 14 FR | Straight | 201-1614 |
| 16 FR | Straight | 201-1616 |
| 10 FR | Coude | 202-1610 |
| 12 FR | Coude | 202-1612 |
| 14 FR | Coude | 202-1614 |
| 16 FR | Coude | 202-1616 |



CompactCath® Lite

Same High Quality At An Affordable Price Point

100% DEHP, BPA, and Latex Free

| Size | Tip | Item Number |
|-------|----------|-------------|
| 12 FR | Straight | 1011-1612 |
| 14 FR | Straight | 1011-1614 |
| 16 FR | Straight | 1011-1616 |
| 14 FR | Coude | 1012-1614 |
| 16 FR | Coude | 1012-1616 |



Toll-Free: 1-888-933-2284 (WEE-CATH), ext 3

www.compactcath.com

CompactCath®

Now Available in Straight & Coudé

Easy, discreet, and complete catheterization at your fingertips. CompactCath® is the only 16" catheter that is pre-lubricated, 100% non-touch, and fits in the palm of your hand – literally. 100% DEHP, BPA, and Latex Free. Designed to fit seamlessly into your life, providing maximum privacy, portability, and convenience.

The CompactCath® Advantage

Pre-Lubricated

CompactCath® is pre-lubricated with silicone oil, offering a lubricious non-sticky, no-drip coating. Less friction means less pain and discomfort during catheterization.

Fire-Polished Eyelets

CompactCath®'s super smooth eyelets offer pain free catheterization, minimizing risk of developing micro-abrasions which can lead to buildup of scar tissue in the long term.

Discreet Disposal

CompactCath®'s retraction dial allows the catheter to be coiled back in place for easy, discreet disposal after use and without making a mess.

Drainage Control

Only CompactCath® offers a unique flow control technology, giving you complete control of when and where to void.

Hygienic, Non-Touch Insertion

The advancer sheath and plastic case protect the catheter from being touched, increasing sterility and reducing risk of potential contamination developing UTI's.

New! CompactCath® Coudé Featuring The World's First "Case Up, Tip Up" Design

As long as the case is held in the vertical position, the tip will also be in the proper upright position, eliminating the need for guide stripes, simplifying insertion and allowing coudé users to catheterize correctly and with confidence.

OneCath®

OneCath® is a high quality catheter at an affordable price point, designed to improve the lives of patients and help in an industry with continuously rising healthcare costs.

Features:

- 100% Fire Polished Eyelets
- Flexi-rigidity
- Noise-Free, Easy-Open Packaging
- Uncoated
- 100% DEHP, BPA, and Latex Free
- Available in a variety of sizes and tips
- Conveniently packaged in boxes of 25

CompactCath® LITE

CompactCath® LITE is the world's smallest full-sized catheter and is 100% touch-free due to a protective sleeve that is covering the entire catheter. Like CompactCath, is lubricated with silicone oil that has proven antimicrobial qualities..

Features:

- Protective sleeve enables easy insertion and helps reduce contamination
- Pre-lubricated with Silicone oil that is known to have antimicrobial qualities
- 100% Fire polished eyelets
- 100% DEHP, BPA, and latex free
- Available in sizes 12FR, 14FR, and 16FR
- Available in Straight and Coude tips

